

URGENT Field Safety Notice (Expansion)

M5071A Adult and M5072A Infant/child pads cartridges for use with HS1/OnSite/Home AEDs may experience gel separation and reduction of gel surface area

15-NOV-2022

This document contains important information for the continued safe and proper use of your equipment

Please review the following information with all members of your staff who need to be aware of the contents of this communication. It is important to understand the implications of this communication.
Please retain a copy with the equipment Instructions for Use.

Dear Customer,

On March 4, 2022, Philips voluntarily initiated an "Urgent Field Safety Notice" for the M5071A Adult and M5072A Infant/child pads cartridges for use with HS1/OnSite/Home AEDs. Due to customer feedback, we have decided to notify customers who received an HS1/OnSite/Home AED more than 10 years ago. **Please respond, even if you no longer own the AED.** This notice has been reported to the appropriate Regulatory Agencies. This letter is intended to inform you about:


1. What the problem is and under what circumstances it can occur

Some electrode pads used with the HS1/OnSite/Home AED have been observed to experience gel separation from the foam/tin backing when peeled from the yellow plastic liner. The gel may fold onto itself, resulting in reduced surface area of gel on the pad. A pad in this condition could cause the HS1/OnSite/Home AED to deliver less effective or ineffective therapy to the patient due to the reduced surface contact area with the skin.



Areas of missing gel

Separated, folded gel may also have a discolored or melted appearance. While the gel may also have a discolored or melted appearance, the appearance does not have any impact on the delivery of therapy; however, there may be a delay in therapy if the user hesitates to apply the pad due to its appearance, and the HS1/OnSite/Home AED may deliver less effective or ineffective therapy to the patient due to the reduced contact area with the skin.



Gel almost completely missing

It is also possible that the gel could separate almost completely from the foam/tin backing when peeled. Due to a small amount of gel surface contact area with the skin, electrical arcing could occur when a shock is delivered leading to burns to the patient's skin, or the HS1/OnSite/Home AED could be unable to deliver any shock through the pads. A delay in therapy will result while the user installs a replacement pads cartridge (if available) or performs CPR while waiting for Emergency Medical Services Personnel to arrive.

2. Hazard/harm associated with the issue: An electrode pad that experiences gel separation could result in less effective therapy for a patient, delay of therapy for a patient, or cause the HS1/OnSite/Home AED to be unable to deliver any shock through the pads.

3. Affected products and how to identify them: M5071A Adult and M5072A Infant/child pads cartridges with a LOT number that begins with "Y" may experience the issue.

LOT Number of affected pads begins with "Y"



4. Actions that should be taken by the customer/user in order to prevent risks for patients or users

Keep your HS1/OnSite/Home AED in service until you receive updated pads. If you need to use your HS1/OnSite/Home AED before an updated pads cartridge has been installed, ensure the majority of the pad surface is covered with gel and apply the pads to the patient. If you notice the gel beginning to separate from the white foam backing as you peel, try if possible to prevent the gel from folding onto itself. Do not hesitate to apply the pads to the patient unless the gel has almost completely separated from the backing. In case of trouble, install a spare pads cartridge if available and continue the rescue. No matter the state of the pads, follow the voice prompts because the HS1/OnSite/Home AED will guide you through the necessary actions. If the problem continues and you do not have a spare pads cartridge, attend to the patient, providing CPR if needed, until Emergency Medical Services Personnel arrive. Philips recommends that you store a spare pads cartridge with your HS1/OnSite/Home AED. A short video showing how to replace the pads cartridge can be found at: www.philips.com/replace-aed-pads

If your HS1/OnSite/Home AED is in service and installed with an unexpired M5071A Adult pads cartridge with a LOT number that begins with "Y", then you are eligible to receive an updated Adult pads cartridge, free-of-charge. Unexpired M5071A spare pads cartridges with a LOT number that begins with "Y" will be replaced, free-of-charge. You must respond to receive any free-of-charge updated Adult pads cartridges. Please pass this notice to all those who need to be aware within your organization or to any organization where HS1/OnSite/Home AED devices or pads cartridges have been transferred. Please keep a copy of this letter with the Instructions for Use/Owner's Manual. **Please respond to this notice because your response is necessary to ensure effectiveness of the recall notification.** Even if you transferred your device to someone, or if your device is no longer in service, please respond.

5. Actions planned by Philips to correct the problem

Philips will provide to customers who respond one updated Adult pads cartridge, free-of-charge, per HS1/OnSite/Home AED in service that is installed with an unexpired M5071A Adult pads cartridge with a LOT number that begins with "Y". Unexpired M5071A spare pads cartridges with a LOT number that begins with "Y" will be replaced, free-of-charge. You must respond to receive any free-of-charge updated Adult pads cartridges.

Infant/child pads cartridge updates will be handled separately. If you own an M5072A Infant/child pads cartridge, Philips will provide, free-of-charge, updated M5072A Infant/child pads cartridges when available to replace unexpired Infant/child pads cartridges.

A similar notice was previously sent to customers who purchased HS1/OnSite/Home AEDs less than 10 years ago. **If you received that notification, please respond to both notifications.** If you need any further information or support concerning this issue, please contact your Philips representative. Philips regrets any inconvenience caused by this problem.

Sincerely,



Tanya DeSchmidt
Director of Quality

Reply Form for FSN C&R 2022-CC-EC-011

Instructions: Whether your HS1/OnSite/Home AED remains in service or not, please respond within the next 30 days. < Add text and/or graphics to suit the method of response you want, e.g., QR code, URL, email, fax, etc. Example: "Please write your answers below, provide your signature and contact info, and return the form. You may take a photo of the completed form or scan with a scanner. Email or fax to: *Insert Philips local Key Market email address and fax number here.*" >

1. What is your name and the name of your facility?

Write your name and the name of your facility.

2. What is the status of your HS1/OnSite/Home AED?

Please check the box and use the space on the right to provide additional info. If different status for multiple AEDs please provide information for all AEDs and make the status of each clear in your response.

My HS1/OnSite/Home AED remains in service, and my M5071A Adult pads cartridge is unexpired and has a LOT number that begins with "Y". The **Serial Number** of my HS1/OnSite/Home AED and the **LOT number(s)** of my M5071A Adult pads cartridge, and spare cartridge(s) if applicable, are listed to the right. (If multiple HS1/OnSite/Home AEDs and Adult pads, please include all. Add a separate page if needed.) Therefore, please send me the replacement Adult pads cartridge(s). I understand I need to provide both the serial number of the device and LOT number(s) of the unexpired Adult pads cartridge(s) to receive the replacement(s) free-of-charge. (Note: M5072A Infant/child pads cartridge updates will be handled separately, so please do not include Infant/child pads cartridge information.)


My HS1/OnSite/Home AED is no longer in service (discarded, lost, permanently retired or destroyed). If the Serial Number is known to me, I have listed it to the right. (If multiple, please include all serial numbers.)

Ownership of my HS1/OnSite/Home AED was transferred to someone else. If available, I have listed the Serial Number of the AED to the right. (If multiple, please include all serial numbers. Add a separate page if needed.) I will share this notice with the new owner(s).


My HS1/OnSite/Home AED is installed with an Adult pads cartridge that is past its expiration date. I have listed the AED Serial Number(s) in the space provided, but I understand that because my pads are expired, I will not receive a free-of-charge replacement Adult pads cartridge.

Use this space to provide additional information.

Serial number example



LOT # example



I acknowledge receipt and understanding of this Urgent Field Safety Notice and confirm that the information from this notification has been properly distributed to all who handle the HS1/OnSite/Home AED. The contact information below will be used to update Philips' records in case of future safety notices, and to provide the replacement pads cartridge(s). The information below will not be used for marketing purposes. Acknowledging this notice will not reduce your coverage or rights under any Philips AED Warranty or Indemnification.

Signature: _____

Printed Name: _____

Address:

(Facility address or personal address, as applicable)

Telephone Number: _____

Email Address: _____

Date (DD-MMM-YYYY): _____