

**MAXIMUM QUANTITY ALLOWED FOR IMPORTATION OF MEDICAL  
DEVICE FOR PERSONAL USE**

<b>NO.</b>	<b>Examples of medical device</b>	<b>Maximum quantity allowed</b>	<b>Remarks/conditions</b>
1	Disposable contact lens	3 months'supply (e.g. 20 pcs/box= 90 pcs)	Depends on usage (monthly/daily)
2	Blood glucose monitoring	1 unit	Household or personal health monitoring
3	Hearing aids	1 unit or 1 pair	Must not be resold or used for clinics
4	Physiotherapy devices	1 unit	Must not be professional or clinical use
5	CPAP	1 unit	Requires medical justification (e.g. sleep apnea diagnosis)
6	BP monitor	1 unit	Household or personal health monitoring
7	Thermometer	1 unit	Household or personal health monitoring
8	In-Vitro Diagnostic (IVD) products	1 piece of each type	
9	In-Vitro maintenance test strips	Quantity for 3 month's supply	

NOTE: The above list is non-exhaustive and may be updated by the Authority from time to time.